



# TEST DIRECTORY



◀ [Back to Test Directory](#)

## TH69 : NOVEL CORONAVIRUS COVID-19 OROPHARYNX

### INFORMATION

**Alternate Name:** COVID-19 Oropharynx

**Clinical Utility:** For the detection of the novel COVID-19, Coronavirus.

**Methodology:** Real-Time Polymerase Chain Reaction

# BioReference

LABORATORIES

an **OPKO** Health Company**Turnaround Time:** 3 days MENU**Preferred Specimen:** Swab-Star Swab**Alternative Specimen:** Swab-Viral Culturette**Specimen Comment:** Swab must be labeled with source: Oropharynx**Storage Instruction:** Refrigerate

## BILLING

**CPT Codes:** 87635 (1)

### Choose Reliability

BioReference is a participating provider in the UnitedHealthcare Preferred Laboratory Network.

[LEARN MORE](#)

TT

**Your feedback is important to us.** If you are you a current customer, please take the time to fill out our Client Satisfaction Survey.

[TAKE SURVEY](#)

**Did you know we are open for visits by physicians and office staff?** Please contact us if you'd like to take a tour of our facilities.

[LAB TOURS](#)



**TEST DIRECTORY**

**VISIT OUR OTHER SITES**



Score®



GeneDx



GenPath Oncology



GenPath Women's Health



Laboratorio Buena Salud



my.labology

**CAREERS**

Visit our careers website for a listing of current openings [here](#).



## STAY CONNECTED



Sign up to get latest updates [here](#).



[EMPLOYEE LOGIN](#)

Copyright 2020 © **BioReference Laboratories, Inc.**

